


# Medix Nebulisers: **AC 4000**



AC 4000 ensure ensures fast and effective therapy everytime for bronchodilators (relievers), corticosteroids (preventers) and antibiotics

- Ideal for home and General Practice use
- Easy to use
- Suitable for patients of all ages
- Discreet storage compartment allows for all accessories and cable to be neatly stored
- Durable and reliable
- Easy carrying handle
- Comes complete with starter pack that contains necessary accessories ready for immediate use (Nebuliser chamber, driveline, inlet filter, adult face mask)
- Shoulder bag available



[www.clement-clarke.com](http://www.clement-clarke.com)



**CLEMENT CLARKE INTERNATIONAL**

Engineered for Life

TOTAL RESPIRATORY SOLUTIONS

Purchase this on our E-Shop at [www.clement-clarke.com/shop](http://www.clement-clarke.com/shop)



## Technical Specifications

Dimensions: 375 x 235 x 115mm

Weight: 2.5kg

Electrical Requirements: 220-240V 50-60Hz

Typical Dynamic Flow: 7 L/min @ 100 Kpa

Particle Size: 3.1 microns\*

% Particles  $\leq$  5 $\mu$ m: 75%\*

Nebulisation Time: 4.5 - 7.5 minutes for 2.5ml

Nebulisation Rate: 0.07ml per minute\*

Compliance: Electrical Safety Class II Type BF Class  
IIa MDD 93/42/EEC CE 0120.

\*Measured in accordance with BSEN 13544-1

## Accessories Included

1 x Microneb III nebuliser chamber, 1 x driveline, 1 x adult face mask, 1 x inlet filter

## Part Codes

AC 4000 (UK Plug): 3604950A

AC 4000 (EU Plug): 3604951

## Consumables



## Year Packs

Clement Clarke advises you change your Nebuliser consumables every 3 months. A Year Pack gives the average user enough consumables for one year.

A year pack contains 4 Microneb III nebuliser chambers, 4 drivelines, 4 adult face masks (or child face masks), 1 mouthpiece, 4 inlet filters and 1 outlet filter.

## Part Codes

Adult Year Pack: L3605122

Child Year Pack: L3605123

Shoulder Bag: 3607691

**Clement Clarke** are well established in offering accurate, precise and dependable products for the diagnosis, monitoring and treatment of respiratory disease

Clement Clarke continues to develop broad product ranges to strengthen our portfolio in **Asthma Management, COPD, Smoking Cessation, Cystic Fibrosis and Allergic Rhinitis** in providing the total respiratory solution.