

Peak Flow: Mini-Wright AFS Low Range



Having set the standard with the Mini-Wright, the AFS Low Range brings the same accuracy and consistency to those patients with lower peak flow monitoring requirements.

- Ideal for children and elderly or those with low Peak Flow measurements
- Incorporates unique Aerodynamic Flow System (AFS)
- Easy to read scale
- Available with windmill trainer
- Compact enough to fit discreetly into bags
- Readings can be taken "on the move"
- Simple to use
- Range of consumables available



www.clement-clarke.com



**CLEMENT CLARKE
INTERNATIONAL**

Engineered for Life

TOTAL RESPIRATORY SOLUTIONS

Purchase this on our E-Shop at www.clement-clarke.com/shop



Technical Specifications

Dimensions: 17.6 x 9.5 x 4 cm

Weight: 96g

Measurement range

30-400 L/min (EU Scale)

30 – 400 L/min (ATS)

Accuracy: +/- 10% or 10 L/min

Repeatability: \pm 5 l/min

Part Codes

Mini Wright AFS Low Range EU Scale: 3104710

Mini Wright AFS Low Range ATS Scale: 3104050

Mini Wright AFS Low Range EU Scale DOH: 3104708

Mini-Wright Windmill Trainer (Pack of 5): 3122015

Consumables



One Way Valve Mouthpieces

Allow Peak Flow Meters to be used as a multi-patient device. Helps to prevent cross infection.

Part Codes

Adult Disposable Expiratory One Way Valve Mouthpieces

Box of 100 (5x20): 3122062

Box of 200: 3122064

Box of 500: 3122066

Paediatric Adaptors are also available for One Way Valve Mouthpieces

Pack of 1: 3122049

Pack of 20: 3122120



Two Way Mouthpieces

Cost effective mouthpieces that allow Peak Flow Meters to be used as a multi-patient device. Cleaning is recommended between each patient when using these mouthpieces.

Part Codes

Adult Disposable Two Way Mouthpieces Box of 100 (5x20): 3125030

Box of 500: 3122003

Paediatric Disposable Two Way Mouthpieces

Box of 40: 3125005

Box of 500: 3123030

Clement Clarke are well established in offering accurate, precise and dependable products for the diagnosis, monitoring and treatment of respiratory disease

Clement Clarke continues to develop broad product ranges to strengthen our portfolio in **Asthma Management, COPD, Smoking Cessation, Cystic Fibrosis and Allergic Rhinitis** in providing the total respiratory solution.

www.clement-clarke.com

DISTRIBUTED BY: