

Spacer Devices: **Able Spacer**



The Able Spacer is a small-volume spacer device which is ideal for people on the move.

- Fits all Metered Dose Inhalers (MDI)
- Assists in drug delivery to the lungs
- Built-in coaching signal maintains correct Inspiratory rate for maximum drug delivery
- Inhaler can be stored inside when not in use
- Ideal for active patients
- Suitable for adults and children
- Easy to use
- Detachable masks available in 3 sizes (Small, Medium and Large)



www.clement-clarke.com



**CLEMENT CLARKE
INTERNATIONAL**

Engineered for Life

TOTAL RESPIRATORY SOLUTIONS

Purchase this on our E-Shop at www.clement-clarke.com/shop



Technical Specifications

Volume: 135ml
 Whistle activation flow:
 50 l/min average +/- 10 l/min
 Resistance to flow:
 Less than 1.6 inches of water @ 30 litres/ second

Part Codes

Able Spacer in a Box: 3607000

Consumables



Able Spacer Masks

Masks are available that are specifically designed for use with the Able Spacer. Comfortable for the patient and easy to connect. These masks are available in Silicon or PVC forms.

Part Codes

- Small Silicon Mask (Infant): 3607005
- Medium Silicon Mask (Child): 3607006
- Large Silicon Mask (Adult): 3607007
- Small PVC Mask (Infant): 3607010
- Medium PVC Mask (Child): 3607009
- Large PVC Mask (Adult): 3607031

Clement Clarke are well established in offering accurate, precise and dependable products for the diagnosis, monitoring and treatment of respiratory disease

Clement Clarke continues to develop broad product ranges to strengthen our portfolio in **Asthma Management, COPD, Smoking Cessation, Cystic Fibrosis and Allergic Rhinitis** in providing the total respiratory solution.

www.clement-clarke.com

DISTRIBUTED BY: