

# AirSep® FreeStyle® 5

PORTABLE OXYGEN CONCENTRATOR

## FREEDOM HAS A NEW LOOK

CAIRE introduces important updates to its high capacity portable oxygen concentrator – the AirSep® FreeStyle® 5. The new design adds features and benefits that will not only make the device more appealing to patients, but even more reliable for providers.

- **New Power Supply Connection**

Now featuring multi-directional installation with no pins to align and no release button to press, it is easier than ever to re-charge on the go.

Robust metal connection points ensure the highest durability of both the power supply and FreeStyle 5 for patients and provider alike.

- **New User-Replaceable Battery**

Swap battery packs on the go instead of waiting to re-charge an internal battery or having to add heavy external batteries. The battery pack will re-charge anytime AC or DC power is connected and is easily accessible from the front of the device.

- **Longer Battery Life**

The new battery pack provides an additional 30 minutes of duration at setting 2, providing up to 2 ½ hours with a single battery pack at setting 2.

- **New Color Scheme**

New color scheme is more appealing and non-medical in appearance. The matte finish will help providers by resisting everyday scratches.



  
**CAIRE**<sup>®</sup>  
A Chart Industries Company

770.721.7759 | Toll Free 1.800.482.2473

[www.CAIREMedical.com](http://www.CAIREMedical.com)



\*\* SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.

PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.

© 2016 CAIRE Inc. All Rights Reserved ML-CONC0079 Rev B

## PROVIDER FEATURES

- Settings 1-5 covers more patients as their disease progresses.
- Non-delivery modality eliminates expensive cylinder deliveries.
- Time-proven and trusted quality
- Standard 3-year warranty includes sieve beds

## PATIENT FEATURES

- User-replaceable battery pack with up to 2½ hour battery life at setting 2
- Easily wearable with a handle, shoulder strap, or backpack
- FAA Approved
- Fully functional on AC or DC power.
- Simple, One-Touch Operation

## CLINICAL FEATURES

- UltraSense™ Technology offers sensitive breath detection
- Quickly detects and pulses to promote oxygen delivery in the earliest part of the breath
- Adjusts bolus size based on breath rate
- Up to 1000 mL/min oxygen output



### SPECIFICATIONS

<b>Flow Rates</b>	Pulse 1-5
<b>Dimensions</b>	10.7 x 6.6 x 4.4 in (27.2 x 16.8 x 11.2 cm)
<b>Weight (w/Battery)</b>	6.7 lbs (3.0 kg)
<b>AC Power</b>	100-240 VAC, 50-60Hz (2 amp max)
<b>DC Power</b>	12-24 VDC (15 amp max)
<b>Operating Temp.</b>	41°F to 104°F (5°C to 40°C)
<b>Sound Level</b>	48 dB(A) @ Setting 3
<b>Operating Altitude</b>	Up to 12,000 ft (3657.6 m)
<b>O2 Concentration</b>	90% (+5.5% / -3%)
<b>Battery Duration</b>	Setting 1 - 4.25 hours (4:15) Setting 2 - 2.5 hours (2:30) Setting 3 - 2.0 hours (2:00) Setting 4 - 1.75 hours (1:45) Setting 5 - 1.25 hours (1:15)
<b>Battery Re-Charge</b>	4 - 5 hours
<b>Warm-Up Time</b>	< 2 minutes
<b>O2 Output</b>	Setting 1 - 227.5 mL/min Setting 2 - 420 mL/min Setting 3 - 630 mL/min Setting 4 - 840 mL/min Setting 5 - 1000 mL/min
<b>Limited Warranty</b>	3 years unit ** 1 year battery & power supplies

### ORDERING INFORMATION

**FreeStyle 5 Includes:** Portable Concentrator, Battery Pack, AC Power Supply, AC Power Cord, DC Power Supply, Carrying Case, Shoulder Strap, User Manual, Carry-all Accessory Bag

<b>AS077-101</b>	USA
<b>AS077-102</b>	European
<b>AS077-B101</b>	USA - 2 batteries
<b>AS077-B102</b>	Euro - 2 batteries
<b>AS077-103</b>	Australia
<b>AS077-104</b>	UK
<b>AS077-B104</b>	UK - 2 batteries
<b>AS077-105</b>	Brazil
<b>AS077-106</b>	Swiss

*\*other countries' power supplies may be available - contact customer service\**

### ACCESSORIES

<b>BT033-1</b>	Spare Battery Pack
<b>MI284-1</b>	Backpack Harness
<b>MI420-1</b>	Carrying Case
<b>MI372-1</b>	Carry-All Bag
<b>MI426-1</b>	Shoulder Strap
<b>F1194-1</b>	Intake Filter
<b>PW034-1</b>	DC Power Supply
<b>PW035-1</b>	AC Power Supply (US)
<b>PW035-2</b>	AC Power Supply (EU)
<b>PW035-3</b>	AC Power Supply (UK)

*Cords Included*

<b>BT032-1</b>	Desktop Charger (US)
<b>BT032-2</b>	Desktop Charger (EU)
<b>BT032-4</b>	Desktop Charger (UK)

*Power Supply Included*

**CAIRE**<sup>®</sup>  
A Chart Industries Company

770.721.7759 | Toll Free 1.800.482.2473

[www.CAIREMedical.com](http://www.CAIREMedical.com)



**\*\* SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.**

**PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.**

© 2016 CAIRE Inc. All Rights Reserved ML-CONC0079 Rev B