

Medix Nebulisers: Econoneb



The Econoneb is a powerful compressor that is continuously rated, making it a cool running nebuliser for demanding use. Its lightweight design makes it ideal for portable use in a hospital.

- Designed for continuous use
- Remains cool during constant use
- Suitable for Hospital, clinic and home use
- Ensures fast and effective therapy every time for bronchodilators (relievers), corticosteroids (preventers) and antibiotics
- Powerful compressor
- Easy to use
- Easy carrying handle
- Reliable
- Comes complete with starter pack that contains necessary accessories ready for immediate use (Nebuliser chamber, driveline, inlet filter, adult face mask)



www.clement-clarke.com



**CLEMENT CLARKE
INTERNATIONAL**

Engineered for Life

TOTAL RESPIRATORY SOLUTIONS

Purchase this on our E-Shop at www.clement-clarke.com/shop



Technical Specifications

Dimensions: 210 x 185 x 185mm

Weight: 2.35kg

Electrical Requirement: 220-240V 50 Hz (110-240V
60 Hz AC model also available)

Typical Dynamic Flow: 8 L/min @ 137 Kpa

Particle Size: 3.1 microns*

% Particles: \leq 5 μ m: 72%*

Nebulisation Time: 2 - 5 minutes for 2.5ml

Nebulisation Rate: 0.01ml per minute*

Compliance: Electrical Safety Class II Type BF Class
IIa MDD 93/42/EEC CE0120

* Measured in accordance with BSEN 13544-1

Accessories Included:

1 x Microneb III nebuliser chamber, 1 x driveline,
1 x adult face mask, 1 x inlet filter

Part Codes

Econoneb (UK Plug): 3605040HW

Econoneb (EU Plug): 3605042HWEU

Consumables



Year Packs

Clement Clarke advises you change your Nebuliser consumables every 3 months. A Year Pack gives the average user enough consumables for one year. A year pack contains 4 Microneb III nebuliser chambers, 4 drivelines, 4 adult face masks (or child face masks), 1 mouthpiece, 4 inlet filters and 1 outlet filter.

Part Codes

Adult Year Pack: L3605122

Child Year Pack: L3605123

Clement Clarke are well established in offering accurate, precise and dependable products for the diagnosis, monitoring and treatment of respiratory disease

Clement Clarke continues to develop broad product ranges to strengthen our portfolio in **Asthma Management, COPD, Smoking Cessation, Cystic Fibrosis and Allergic Rhinitis** in providing the total respiratory solution.

DISTRIBUTED BY: