

NIDEK *Medical*

Proudly Presents
The Newest and Most
Innovative Member
of the Nuvo Family
of Oxygen Concentrators



.....Mark 5 **Nuvo** *Lite*

- "Lite" - a mere 32 lbs. (14.5kg)
- Quiet Operation
- Compact
- Lockable Flow Control Valve
- Liter Flow -- 0.125 - 5 LPM



Metal Outlet Fitting



Lockable Flow Control Valve



Resettable Circuit Breaker

Models

Nuvo Lite Oxygen Concentrator

115 Volt - 60 Hz

520

230 Volt - 50 Hz

920

Nuvo Lite Oxygen Concentrator with Oxygen Concentration Status Indicator (OCSI)

525

925

Specifications*

Average Power Consumption
Oxygen Concentration (at 5 LPM)
Liter Flow
Sound Level
Outlet Pressure
Weight
Dimensions

330 Watts
90% (+6.5% / -3%)
0.125 to 5 LPM
40 dBA typical
7.0 PSI
32 lbs
14"W x 9"D x 23"H

300 Watts
90% (+6.5% / -3%)
0.125 to 5 LPM
40 dBA typical
7.0 PSI
14.5 kg
36 cm x 23 cm x 58.5 cm

Accessories and Re-order Items

Compressor Intake Filter
External Cabinet Filter
Humidifier
Cannula with 7' tubing
Nuvo Lite User's Guide

8400-1180
8400-1025
9012-8774
9012-8780
2010-8400

8400-1180
8400-1025
9012-8774
9012-8780
2010-8401CE (OCSI)

* NOM - Specifications are subject to change without notice

2010-8400LS Rev B