

FreeStyle™ 5

NEW!

Life to the Fullest – Higher Capacity POC

FreeStyle™ 5* is a new, full-capacity portable oxygen concentrator (POC) from AirSep – the world's expert in oxygen generating technology.

The FreeStyle 5 POC offers 5 pulse setting options, which ensure that even your more challenging ambulatory O₂ patients can exercise and have access to the oxygen they need – all in a lightweight, compact, portable unit.

FreeStyle 5 is an attractive, contemporary, non medical-looking device whose pulse-flow delivery enables higher-flow patients to be as active as possible. It's the lightest POC in its class.

Easy to Use, Easy to Operate

The FreeStyle 5 concentrator is easy to use and easy to operate. This model provides patients with all the advantages of the standard FreeStyle unit, including proven reliability and durability in a slightly larger platform for greater capacity. In fact, FreeStyle 5 is the lightest higher-capacity POC on the market at 5.8 lbs (2.6 kg).

Multiple Power Options

FreeStyle 5 can be conveniently used and recharged simultaneously through a universal AC/DC power supply. This allows a patient to depend on the flexibility of what FreeStyle 5 offers, whenever and wherever the oxygen unit is used.

A lightweight AirBelt™ accessory, which was designed to be conveniently worn around the waist, ensures that extended battery time between recharges is possible. This option can be marketed to patients, affording them the ability to maximize duration time on FreeStyle 5 for those occasions when AC or DC power is not readily available.

The POC keypad and its AirBelt both show a digital display of total remaining battery charge in 25 percent increments, which corresponds to how much operating time remains before needing to recharge the unit.



**Digital Keypad/
Alarm Display**

Full 5 LPM equivalency for FreeStyle 5



SPECIFICATIONS

Oxygen Concentration**:	1-5 pulse settings; equivalent to a continuous flow of 1-5 LPM
Dimensions:	10.5 in. H x 6.4 in. W x 4.4 in. D (26.7 cm H x 16.3 cm W x 11.2 cm D)
Weight:	5.8 lb (2.6 kg)
Power:	100-250 VAC, 50/60 Hz; 11-16 VDC; rechargeable, lithium battery
Battery Duration:	Up to 2 hours on internal battery; up to 7.5 hours total when used with lightweight AirBelt™
Audible Alarms and Visual Indicators:	Low battery Cannula disconnect System overdraw High and low pressure High temperature

**Based on an atmospheric pressure of 14.7 psia (101 kPa) at 70° F (21 °C).

Note: Specifications subject to change without notice.

INTERNATIONAL APPROVALS

AirSep Corporation is ISO 9001:2008-certified.

ORDERING INFORMATION

Contact your AirSep representative for availability.



AirSep Corporation
Medical Products Division
401 Creekside Drive
Buffalo, NY 14228-2085 USA

Tel: 716.691.0202
Fax: 716.691.4141
E-mail: mpd@airsep.com
Internet: www.airsep.com

*Product not for sale. For demonstration purposes only.