

# WYMEDICAL

S U R G I C A L S T O C K I N G S

WYMEDICAL

WYMEDICAL

S U R G I C A L S T O C K I N G S


S U R G I C A L S T O C K I N G S

Below knee and open toe surgical  
weight hose for men and women

Graduated compression and  
Physician recommended.

Clinically proven relief  
for tired, aching legs,  
varicosities and edema.

Recommended wear  
during pregnancy.



WYMEDICAL

S U R G I C A L  
S T O C K I N G S

When your physician recommends an under knee style of support stocking with firm compression, choose **WyMedical** Surgical Support Stockings . This style of stocking can help relieve problems associated with:

**Leg Pains** - such as cramps or aches caused by restricted circulation associated with varicose veins.

**Pregnancy** - when extra pressure in the lower abdomen restricts blood flow, requiring assistance in maintaining circulation.

**Post Surgery** - when firm support stockings are recommended for leg veins and muscles weakened by surgery or lack of exercise.

**Professional/Occupational Hazard** - restricted circulation caused by long periods of sitting; sluggish blood flow caused by hours of standing resulting in tired, aching legs, or varicose veins.

Support hosiery, to be effective, must be designed with medically accepted graduated compression.. This compression, strongest at the ankle, gradually reducing towards the top of the stocking, assists blood circulation up the leg and helps to relieve leg discomfort due to swelling, varicose veins, fatigue, pregnancy and other related problems.

How to select your correct size								
Length, floor to knee	Calf Measurement							
	28 cm	30.5 cm	33 cm	35.5 cm	38 cm	40.5 cm	43 cm	45.5 cm
35.5 cm	S	S	S	S	M	M	L	L
38 cm	S	S	S	M	M	L	L	L
40.5 cm	S	S	M	M	M	L	L	L
43 cm	S	M	M	M	M	L	L	XL
45.5 cm	M	M	M	M	L	L	XL	XL
48 cm	M	M	M	M	L	XL	XL	XL

Distributed in Australia by Wymedical (08) 362 1852  
Made in USA