

# P2-E Portable Oxygen Concentrator

Lightweight power. Reliable oxygen.  
Everyday freedom.

The Kingon P2-E is the world's smallest and lightest portable oxygen concentrator. Weighing just 1.98 kg and measuring 22.1 x 8.5 x 16 cm, it's the closest you can get to a pocket-sized POC. Take it anywhere—at home or on the go—and lift it easily with one hand.




-  **Lightweight 1.98 kg** - compact and portable
-  **More than 5 hours battery life** - stay powered longer
-  **Smart pulse-dose delivery** - oxygen only when you inhale
-  **Quiet 37 dB operation** - low noise for discreet use
-  **Airline approved** - travel confidently

## At a Glance

Specification	Details
Oxygen Concentration	90% (-3% / +6%) at all settings
Flow Settings	1-5 (pulse dose delivery)
Weight	1.98 kg (with battery)
Size	22.1 L x 16 H x 8.5 W cm
Noise Level	37 dB (setting 2)
Battery Duration	1: 5h24 2: 4h30 3: 3h18 4: 2h30 5: 2h06
Recharge Time	Approx. 4 hours (AC/DC)
Power Options	AC 100-240 V · DC 12 V · Battery
Display	2.8" LCD screen
Warranty	3 yrs device · 12 mths sieve bed/batteries

## Package Inclusions



Discover the Kingon difference.

Talk to your oxygen provider or visit [kingon.com.au](http://kingon.com.au) to learn more.



## Live Freely. Breathe Confidently.

The Kingon P2-E gives you freedom of movement without the restrictions of heavy oxygen cylinders or large at-home concentrators.

The Kingon P2-E gives you freedom of movement without the restrictions of heavy oxygen cylinders or large at-home concentrators. Whether relaxing at home, heading out for daily activities, or travelling, the P2-E provides reliable oxygen therapy in a highly compact form.

### Simple, User-Friendly Design

The Kingon P2-E features an easy-to-navigate interface with a bright 2.8" LCD screen, ideal for users with impaired vision. Adjustable flow settings allow the device to adapt as your oxygen needs change. Visual alerts remain on-screen even when silenced, ensuring important information—like low battery, low oxygen purity, or no breath detected—is always visible.



#### Reliable Performance

Ideal for home and on-the-go oxygen therapy.



#### Simple Maintenance

Easy-to-clean filters and accessible parts.



#### Quiet Operation

Only 37 dB — as quiet as library ambience.



#### Trusted Support

Comprehensive Australian warranty and service.

Experience uncompromised freedom with Kingon P2-E.

Visit [kingon.com.au](http://kingon.com.au) or speak to your oxygen provider today.