



Live Freely. Breathe Confidently.

The Kingon P2-E6 gives you freedom of movement without the restrictions of heavy oxygen cylinders or large at-home concentrators.

Whether relaxing at home, heading out for daily activities, or travelling, the P2-E6 provides reliable oxygen therapy in a highly compact form.

Simple, User-Friendly Design

The Kingon P2-E6 features an easy-to-navigate interface with a bright 2.8" LCD screen, ideal for users with impaired vision. Adjustable flow settings allow the device to adapt as your oxygen needs change. Visual alerts remain on-screen even when silenced, ensuring important information, like low battery, low oxygen purity, or no breath detected, is always visible.



Reliable Performance

Ideal for home and on-the-go oxygen therapy.



Simple Maintenance

Easy-to-clean filters and accessible parts.



Quiet Operation

Only 37 dB as quiet as library ambience.



Trusted Support






Comprehensive Australian warranty and service.

P2-E6 Portable Oxygen Concentrator

**Lightweight power. Reliable oxygen.
Everyday freedom.**

The Kingon P2-E6 is the world's smallest and lightest portable oxygen concentrator in its class. Weighing just 1.98 kg and measuring 22.1 × 8.5 × 16 cm.. Take it anywhere, at home or on the go, and lift it easily with one hand.



-  **Lightweight 1.98 kg** – compact and portable
-  **More than 5 hours battery life** – stay powered longer
-  **Smart pulse-dose delivery** – oxygen only when you inhale
-  **Quiet 37 dB operation** – low noise for discreet use
-  **Airline approved** – travel confidently

At a Glance

Specification	Details
Oxygen Concentration	90% (-3% / +6%) at all settings
Flow Settings	1 - 6 (pulse dose delivery)
Weight	1.98 kg (with battery)
Size	22.1 L × 16 H × 8.5 W cm
Noise Level	37 dB (setting 2)
Battery Duration	Up to 5 hours 24 mins duration on setting 1
Recharge Time	Approx. 4 hours (AC/DC)
Power Options	AC 100–240 V · DC 12 V · Battery
Display	2.8" LCD screen
Warranty	3 yrs device · 12 mths sieve bed/batteries

Package Inclusions



Distributor:
WyMedical
 466 Fullarton Road, Myrtle Bank, SA, 5064
 Toll Free: 1800 812 097
 E: sales@wymedical.com.au W: www.wymedical.com.au